
What is included with a standard TOOsonix system?

The research systems from TOOsonix consist of the base structure to enable HIFU treatments directly after delivery and training of relevant users. The system includes the following units and functions:

- UPU - Ultrasonic Power Unit.
 - Frequency range 5 - 30 MHz
 - Universal power input
- 2 Pcs Handpieces
 - TOOsonix HPU-001 handpiece
 - 2 m cable with composite connector
 - Resonance Frequency 20 MHz
 - Nominal Focal Dept 1.7 mm and 2.7 mm for treatment of upper and lower dermis
 - 2 MP optical vision system
 - 30 kJ Energy credit"
- Medical 22" medical panel PC
 - Medical keyboard (for service use only)
 - Medical Mouse (for service use only)
- TOOsonix Software
 - Graphical User interface for ONE-R
- Mobile Cart
 - Mobile cart incl shelve for storage of consumables etc.
- Footswitch
 - Medical Footswitch for activation of ultrasound
- Ultrasound Tester
 - Oil-water mixture for test of HIFU functionality
- Ultrasound Gel
 - 2 x 250 ml ultrasound gel for acoustic coupling
- DI Water for transducer filling
 - 2 x 500 ml coupling agent for transducer chamber
- Pipette for water transducer filling
 - 10 x Specialized pipette for filling of transducer chamber
- Acoustic Window
 - LDPE film - 30m transparent acoustic window
- O-rings for Acoustic Window
 - 10 x O-rings for fixating the acoustic window in the transducer head

In addition to the above, additional handpieces, for example to cover other frequencies or focal depts can often be relevant. Such handpieces can be ordered separately according to customer specification.

We furthermore welcome that any other special requests are discussed with our engineering department.

Contact

TOOsonix A/S
Agern Allé 1
DK-2970 Horsholm
Denmark

info@toosonix.com
www.toosonix.com

What is included with a standard TOOsonix system?

Frequently Asked Questions

The company policy is that customers should only be satisfied if they have been fully and effectively supported through the entire process from the first inquiry to the final handover of a system to trained and qualified users.

A satisfied customer that has been enabled with an operational system according to their specifications will in turn enable TOOsonix in obtaining our overall mission:

Helping people by use of high frequency HIFU




Fig. 1. Standard TOOsonix HIFU System.